# Cyanoborohydride (BH<sub>3</sub>CN)

Cyanoborohydride Resin, MP, is a macroporous resin functionalized with a tetraalkylammonium cyanoborohydride end group. It is a commonly utilized reagent for the conversion of carbonyl containing moieties to amines via reductive amination protocols. It is also known to reduce a variety of imminuim ions including pyridinium salts. aldehydes and  $\alpha,\beta$ -unsaturated carbonyls to the corresponding alcohols has also been reported.

A standard protocol involves combining requisite carbonyl and amine in a slightly acidic solution, followed by agitation in the presence of MP-Cyanoborohydride. Simple filtration provides a mixture of desired product (as the acid salt) and either excess amine or carbonyl. Scavenging the excess reagent with the appropriate scavenger, followed by free-basing with MP-carbonate, provides only the desired product upon concentration.

Safety concerns arising from cyanide contamination are alleviated since the toxic species is retained on the polymer support for effective containment and disposal.



$$N$$
  $BH_3(CN)$ 

#### General Reaction

$$\begin{bmatrix} H_{\bullet} & R_2 \\ H_{\bullet} & N \end{bmatrix} AcO^{\ominus} \xrightarrow{N} (CO_3^{2-})_{0.5} \qquad H_{\bullet} & R_2$$

#### References

Hutchins, R. O. J. Chem. Soc., Chem. Commun. 1978, 1088-1089. Habermann, J. J. Chem. Soc., Perkin Trans. 1 1998, 3127-3130.

Ley, S. V. J. Chem. Soc., Perkin Trans. 1 1998, 2239-2241.

Chesney, A. Green Chem. 1999, 1, 209-219.

Baxendale, I. R. Bioorg. Med. Chem. Lett. 2000, 10, 1983-1986.

Hemming, K. Synlett 2000, 1565-1568.

Kirschning, A. J. Prakt. Chem. 2000, 342, 508-511.

Blackburn, L. Org. Lett. 2001, 3, 1637-1639.

Contour-Galcera, M. O.; Bioorg. Med. Chem. Lett. 2001, 11, 741-745.

### Solvent Compatibility

THF

**DCM** 

DCE

**DMF** 

**NMP** MeOH

**EtOH** 

## Ordering Information

### MP-Cyanoborohydride

Loading: 2.8-3.0 mmol/g	10g	SPMP 05-10
	25g	SPMP 05-25
Bead size: 330-1225 microns, 15-50 mesh	100g	SPMP 05-100
(>90% within)	1Kg	SPMP 05-1kg

For additional information contact info@suprasciences.com or visit www.suprasciences.com